

Fresas de Diamante

Roentgen®

ISO NO.	246/007	247/007	108/009	171/013	196/020	196/105	197/012	289/013	150/012	237/020	298/013	298/015	254/015	160/015
Model NO.	CD-56F	CD-57F	CD-58F	CD-59F	CR-11F	CR-12F	CR-21F	CR-22F	EX-18F	EX-21F	FO-21F	FO-22F	FO-30F	TC-11F
Head length	2.1	3.5	3.1	6.8	5.0	4.0	7.0	8.0	0.1	5.0	7.0	7.0	3.2	10.0
Total length	16.7	16.7	16.2	16.7	19.0	19.0	21.5	21.5	21.0	18.9	21.0	21.0	19.0	21.9
ISO NO.	160/013	173/015	171/015	199/015	198/015	197/015	199/015	199/017	304/015	304/015	001/016	196/019	196/014	237/019
Model NO.	TC-21F	TF-12F	TF-21F	TR-11F	TR-13F	TR-21F	TR-25F	TR-26F	WR-31F	WR-S31F	BR-40EF	CR-11EF	CR-12EF	EX-21EF
Head length	8.0	10.0	4.0	10.0	9.0	7.0	10.0	9.0	0.55	0.55	1.5	5.0	4.0	5.0
Total length	21.9	21.7	19.7	21.7	21.7	21.0	21.6	21.7	19.1	16.6	19.0	19.0	19.0	18.9
ISO NO.	248/013	298/012	298/014	243/010	248/010	248/011	160/014	160/012	173/014	171/014	199/014	198/015	197/013	199/014
Model NO.	FO-20EF	FO-21EF	FO-22EF	FO-40EF	FO-41EF	FO-42EF	TC-11EF	TC-21EF	TF-12EF	TF-21EF	TR-11EF	TR-13EF	TR-21EF	TR-25EF
Head length	6.1	7.0	7.0	4.3	6.0	7.7	10.0	8.0	10.0	7.0	10.0	9.0	7.0	10.0
Total length	20.0	21.0	21.0	23.5	18.5	22.0	21.9	21.9	21.7	20.9	21.7	21.7	21.0	21.6
ISO NO.	199/016	001/019	257/019	190/019	111/015	001/015	160/017	173/019	198/019	197/018	068/043	655/018	068/042	166/011
Model NO.	TR-26EF	BR-31C	FO-32C	FO-54C	SF-12C	SI-47C	TC-11C	TF-13C	TR-13C	TR-62C	WR-13C	EX-28	WR-11	TC-10
Head length	9.0	1.6	4.0	9.5	8.1	1.7	10.1	10.1	9.1	7.1	1.9	8.0	1.0	10.0
Total length	21.7	19.1	19.7	21.0	22.0	19.0	22.0	21.8	21.8	20.0	19.1	21.5	19.0	22.0


Diamond - Burs series

Notes: every diamond bur's rotation speed should be slower than the allowed max speed according to the above color remarks.

450,000 R.P.M. 300,000 R.P.M. 160,000 R.P.M.

1 name:

- * TF - taper flat
- * SR - straight round
- * WR - wheel round
- * S - Short handle
- * TR - taper round
- * EX - special shape
- * CR - inlay prepared diamond-burs
- * SS - extreme short handle; C-coarse
- * TC - taper cone
- * BR - ball round
- * FO - flame point
- * BC - ball convexity
- * CD - children dental diamond burs
- * F - fine
- * SF - straight flat
- * DI - double inverted cone
- * EF - extreme fine
- * SO - straight point
- * SI - single inverted cone
- * RS - round sector



2 color code:

Green Coarse (C) / 125-150u

Blue Standard 106-125u

Red Fine (F) 53-63u

Yellow Extreme fine (EF) 20-30u